

## Hotgen® Antigen Self Test

**Officially recognized in the EU as a lay test (with CE seal).**

This rapid, single-use test is used for the determination of coronavirus antigen. The sample is taken from a human nasal swab. The Corona Rapid Test is used as a rapid test for suspected coronavirus cases. This test can also be used for reconfirmation for nucleic acid detection in discharged patients. This test has CE approval for lay use (self-testing) for the detection of SARS-CoV-2 and is also listed with the BfArM, as well as evaluated by the Paul Ehrlich Institute.

### Features:

- Fast result - 15-20 minutes
- Easy to use even by laypersons - all in one kit
- Hygienically individually packaged test kits
- High accuracy and variant tested (including Omicron)
- Test with CE certification as a lay test

### Scope of delivery:

- Coronavirus antigen test cassette
- Sample extraction buffer
- Sterile swab
- Disposable pressure bag for disposal
- Instructions for use

**Manufacturer:** Beijing Hotgen Biotech Co, Ltd.

**Instrument:** Can be performed without any other equipment

**Format:** Anterior nasal test

**Storage temperature:** 4 - 30°C

**Selection time:** 15 minutes

**Specificity:** 99,76 %

**Sensitivity:** 97,10%

**CE Certification:** CE 0123

**BfArM-AT number:** AT1236/21

**Packaging:** individually packaged



**CE-Certificate**



**Variant Verification (Omicron)**



**Watch video**

